Therabody

# RecoveryPulse<sup>™</sup> Arm

Vibration Device

Powered by **QX-Micro Motor** Technology<sup>™</sup>

## Table of contents

EN	3 - 9
DE	10 - 16
FR	17 - 23
IT	24 - 30
ES	31 - 37
NL	38 - 44
PT	45 - 51
RU	52 - 58

## Product Overview

# RecoveryPulse Arm

## Vibration Device

*Compression with vibration therapy.*

From the creators of Theragun® comes RecoveryPulse, compression and vibration therapy in one wearable device. From performance prep to post-workout relief to overall wellness, there is something for everyone. Proprietary lightweight motors along the arm enhance warm-up and improve recovery.

## Intended Use


The RecoveryPulse Arm Vibration Device is intended to help reduce soreness and enhance mobility, focusing on the biceps, triceps, and forearm. This wearable vibration device is also intended to increase blood flow to help reduce symptoms commonly associated with acute muscle pain.

## Getting started

### Step 1

#### Put on your wearable device

1. Find the tag with the circular Therabody logo in your arm sleeve. It should land in the center of the shoulder facing outwards once the sleeve is properly placed (See image).
2. Slide your left or right arm into the sleeve, pulling up steadily and gently.
3. Make sure the sleeve is placed correctly: circular Therabody logo tag is centered on your shoulder facing outwards and the battery pocket is easily accessible below the tag. The Therabody logo should be on top of your forearm, next to the outside part of your wrist.
4. Insert the battery in the pocket. Make sure the golden connectors are facing down and the Therabody logo on the battery faces outwards.

5. Firmly press the button below the battery for 2 seconds to turn on the vibration therapy.
6. Press the button again to cycle through our different vibration therapy modes.
7. Firmly press the button again for 2 seconds to turn off the device.

## Step 2

### Feel the Difference.® Every day.

Whether used between workouts or for general relief, take a deep breath and treat yourself as the combination of material, integrated motors, and compression help keep you feeling your best.


See Warnings for important safety instructions

Scan the QR code to get started.

## Getting to know your wearable device

- A** *Button Pad*
- B** *Lithium-ion battery*
- C** *USB-C Port*
- D** *Proprietary lightweight integrated motors*

## Using the RecoveryPulse Arm Vibration Device

1. Put the RecoveryPulse Arm Vibration Device on either the right or left arm. You will know the sleeve is on correctly when the logo is facing up and out. (See Figure 1).
2. Once in a comfortable position, turn ON the integrated vibration motors by firmly pressing and holding the button for 2 seconds.
3. Press the button again to toggle through the range of vibration frequencies. Adjust the vibration frequency as needed.
4. To turn OFF the RecoveryPulse Arm Vibration Device, firmly press and hold the button for 2 seconds.


Figure 1

## Modes Overview

*For general recovery and pain relief*

### Mode 1 - Low:

Vibration therapy at low frequency

*For warm ups and exercise recovery*

### Mode 2 - High:

Vibration therapy at high frequency

*For relaxation and wellness*

### Mode 3 - Wave:

A tailored pattern from low to high frequency to maximize the contrast and stimulation to the nervous system

## Cleaning Instructions:

1. Remove the battery before washing.
2. Place the RecoveryPulse sleeve in the provided laundry bag and machine wash on delicate mode with temperatures below 30°C/80°F.
3. After washing, lay flat to dry. Do NOT machine dry the RecoveryPulse sleeve. Note: Make sure the RecoveryPulse sleeve is completely dry before re-inserting the battery.

Note: Make sure the RecoveryPulse sleeve is completely dry before re-inserting the battery.

## Charging

- The RecoveryPulse Arm Vibration Device is USB-C enabled.
- The RecoveryPulse Arm Vibration Device includes a Therabody USB-C to USB-A cable.
- You can plug the device into any standard USB adapter.
- The USB port is located in the right side of the battery pack and is accessible by removing the battery from the sleeve pocket.
- The battery level is indicated in the LED indicator of the battery:
  - Red: low charge, blue: mid charge, green: fully charged.
  - When charging, the LED will flash to indicate it is properly connected and will turn solid green when complete.
  - You can always check the battery level by inserting the battery in the sleeve and tapping the power button once.

Note: Make sure the charger is from a certified manufacturer and it has not suffered any structural damage. The device will not work while charging.

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